

## BT PULL

BRACE FOR NEUROLOGICAL DISORDERS PARALYSIS OF THE LIMBS WITH HYPERFLEXION OF THE PHALANGES AFFECTING THE SCIATIC NERVE, THE COMMON FIBULAR NERVE, BRACHIAL PLEXUS, OR RADIAL NERVE.



The brace can be used for neurological disorders that cause severe proprioceptive deficits, which mean that the dog is unable to sense the position of the limb on the ground, thus tending to stand on the back of its paw (dorsiflexion). BT Pull is also used in cases of damage to the sciatic nerve (for example as a result of a fracture of the pelvis or femur), in order to ensure that the neurological problem does not lead to orthopaedic damage due to shortening of the flexor tendons of the foot. The brace can be used to: prevent incorrect application of weight on the limb; ensure correct, continuous stimulation of the proprioceptive nerves, also when walking; provide constant stretching of the toes and thus of the flexor tendons of the toes, thereby ensuring that they do not shorten.

### FEATURES

- ▶ Breathable fabric.
- ▶ Hand washable.
- ▶ Excellent wearability.
- ▶ The package contains a reinforcement made of thick material for use in cases of severe hyperflexion of the phalanges, ensuring that they open more completely, also in the case of big dogs. The reinforcement can be removed if not needed.

### RECOMMENDATIONS FOR USE

- ▶ During the period of restraint, it is advisable not to leave the dog alone for too long, because it might try to remove the brace by biting or scratching it.
- ▶ The knot at the top of the brace can be untied so that the elastic fascia can be adjusted to increase or decrease its tension, thus adjusting the opening of the dog's phalanges.

### SIZES

Available in four sizes.



5-10 kg




10-20 kg



20-40 kg



over 40 kg



balto  
ORTHOPAEDIC BRACES FOR ANIMALS